

CPC 250: Purification of [6]-Gingerol from Ginger



TECHNICAL NOTE TN207

G. AUDO AND C. LE QUEMENER | Gilson Purification Applications Laboratory
S. HOLDERITH AND Y. MAESTRO | Chanel Parful Beauté R&T

INTRODUCTION

Ginger, isolated from rhizome of *Zingiber officinale roscoe* (Figure 1) is widely known for its use in flavoring foods and beverages, but it is also commonly used for medicinal, perfumery, and esthetic purposes as well. Natural products from ginger have been shown to exhibit several beneficial bioactivities, including chemopreventative, anti-inflammatory, antioxidant, and anti-emetic.¹ The two main components contributing to the pungent fragrance and flavor exhibited by fresh ginger are gingerols and shogaols, with [6]-gingerol being the major phenol providing ginger's pungent characteristics (Figure 2). [6]-Gingerol is normally obtained as a pungent yellow oil, but can also be found in a crystalline form under the right conditions.

This technical note highlights the ability of a centrifugal partition chromatography (CPC) system (Gilson CPC 250) to obtain high purity fractions of an active compound from a complex crude extract using five times less solvent than flash chromatography and prep HPLC.

MATERIALS AND METHODS

Sample

Three separate preparative runs were performed using 0.5, 1, and 2 grams of crude extract of *Z. officinale roscoe* rhizome.

Apparatus

A Gilson CPC 250 was connected to a Gilson PLC 2050 Purification System equipped with a 50 mL/min quaternary gradient pump, UV/Vis detector, and a fraction collector. The purification runs were performed using Gilson Glider CPC software to control the CPC with the PLC 2050 Purification System (Figure 3).

All HPLC work was performed using the LaChrom Elite HPLC system (VWR) equipped with a photodiode array detector (PDA) (200–800 nm).



Figure 1
Ginger (*Zingiber officinale roscoe*) rhizome




Figure 2
Structure of Gingerol



Figure 3
CPC 250 Connected to a
PLC 2050 Purification System

Analytical HPLC Method

The crude extract was first analyzed by HPLC under the conditions listed in Table 1. [6]-Gingerol is identified at RT=11.17 minutes with a peak area of 19.6% at 210 nm (Figure 4).

Table 1:

Analytical HPLC Method Conditions

HPLC Column	Purosphere RP18 (250*4.6 mm, 5 μm)
Mobile Phase A	Water
Mobile Phase B	Acetonitrile
Gradient	0 min: 40% B 0-55 min: 40% B to 95% B 55-65 min: 95% B 65-70 min: 95% B to 40% B 70-80 min: 40% B
Flow Rate	1 mL/min
Injection Volume	2 μL
Temperature	40°C
Detection	210 nm

CPC Method

The CPC solvent system used was determined by the shake flask method in order to obtain a K_d value near 1.

$$K_d = -1 = \frac{[HPLC \text{ peak area of gingerol}]_{stat}}{[HPLC \text{ peak area of gingerol}]_{mobile}}$$

Once the CPC solvent system was determined, three CPC runs were performed on the Gilson SCPC-250 using 0.5, 1, and 2 g of crude extract. The crude extract sample was dissolved in 5 mL of upper phase and 5 mL of lower phase, filtered through a 0.45 μL membrane filter, and then injected in the CPC according to the conditions described in Table 2. The fractions obtained as a result of the CPC analysis were analyzed by HPLC and grouped according to HPLC purity.




Figure 4
HPLC Analysis of Crude Extract




Figure 5
CPC chromatogram of 500 mg injection




Figure 6
HPLC Analysis of CPC Fraction 15

Table 2

CPC Method Conditions

Column	CPC 250
Column Volume	250 mL
Elution Flow Rate	10 mL/min
Extrusion Flow Rate	30 mL/min
Rotation Speed	2000 rpm
Solvent System	Hept/AcOEt/MeOH/W
Mode	Ascending
Samples	0.5 g crude extract in 5 mL lower + 5 mL upper 1 g crude extract in 5 mL lower + 5 mL upper 2 g crude extract in 5 mL lower + 5 mL upper
Detection	210 nm

RESULTS

ITERATION	EXTRACT INJECTED	RUN TIME	SOLVENT CONSUMPTION	[6]-GINGEROL RECOVERED	YIELD	HPLC PURITY
1	0.5 g	30 min	600 mL	31 mg	6.2%	90%
2	1 g			69 mg	6.8%	96%
3	2 g			87 mg	4.5%	92%

CONCLUSIONS

The combination of the CPC 250 with the PLC 2050 Purification System allows for quick identification and separation of [6]-gingerol from a crude extract. Injections ranging from 0.5–2 g of crude extract were performed to evaluate the change in the amount of [6]-gingerol recovered. The injection of 2 g of crude extract did not significantly decrease the recovery of [6]-gingerol compared to the injection of 0.5 g. CPC provides a rapid one-step method for purification of natural products with minimal solvent consumption.

REFERENCES

1. Khodaie, L. and Sadeghpour, O. (2015) Ginger From Ancient Times to the New Outlook. *Jundishapur Journal of Natural Pharmaceutical Products*, 10(1), e18402.

Trademarks

All product and company names are trademarks™ or registered® trademarks of their respective holders. Use of the trademark(s) in this document does not imply any affiliation with or endorsements by the trademark holder(s).

Notice

This application note has been produced and edited using information that was available at the time of publication. This application note is subject to revision without prior notice.